

New **FDA-Approved** Broad Companion Diagnostic (CDx) for Solid Tumors

## FDA-Approved Content

## Report Section 1

		PATIENT Jane Sample	TUMOR TYPE Lung adenocarcinoma	TRF# TRFXXXXXX
<b>PATIENT</b>	<b>PHYSICIAN</b>	<b>SPECIMEN</b>		
DISEASE Lung adenocarcinoma NAME Not Given DATE OF BIRTH Not Given SEX Female MEDICAL RECORD # Not Given	ORDERING PHYSICIAN Not Given MEDICAL FACILITY Not Given ADDITIONAL RECIPIENT Not Given MEDICAL FACILITY ID Not Given PATHOLOGIST Not Given	SPECIMEN SITE Not Given SPECIMEN ID Not Given SPECIMEN TYPE Not Given DATE OF COLLECTION Not Given SPECIMEN RECEIVED Not Given		
<b>1</b>				
<b>CDx Associated Findings</b>				
<b>GENOMIC FINDINGS DETECTED</b>		<b>FDA-APPROVED THERAPEUTIC OPTIONS</b>		
<b>EGFR</b> L858R		Gilotrif® (Afatinib) Iressa® (Gefitinib) Tarceva® (Erlotinib)		
<b>2</b>				
<b>OTHER ALTERATIONS &amp; BIOMARKERS IDENTIFIED</b>				
Results reported in this section are not prescriptive or conclusive for labeled use of any specific therapeutic product. See professional services section for additional information.				
<b>Microsatellite Status</b> MS-Stable <sup>§</sup>		<b>PTCH1</b> T416S		
<b>Tumor Mutation Burden</b> 11 Muts/Mb <sup>§</sup>		<b>RBM10</b> Q494*		
<b>CDKN2A/B</b> loss <sup>§</sup>		<b>TP53</b> R267P		
<b>EGFR</b> amplification <sup>§</sup>				
<sup>§</sup> Refer to appendix for limitation statements related to detection of any copy number alterations, gene rearrangements, MSI or TMB result in this section.				
Please refer to appendix for Explanation of Clinical Significance Classification and for variants of unknown significance (VUS).				

**1 FDA-Approved Therapies**  
List of FDA-approved companion diagnostics to identify patients who may benefit from associated therapies

**2 All Other Biomarkers**  
All other biomarkers, including tumor mutational burden (TMB) and microsatellite instability (MSI), without companion diagnostic claims

# Professional Services

## Report Section 2



PATIENT  
Jane Sample

TUMOR TYPE  
Lung adenocarcinoma

TRF#  
TRFXXXXXX

---

*Interpretive content on this page and subsequent pages is provided as a professional service, and is not reviewed or approved by the FDA.*

**PATIENT**

DISEASE Lung adenocarcinoma  
NAME Not Given  
DATE OF BIRTH Not Given  
SEX Female  
MEDICAL RECORD # Not Given

**PHYSICIAN**

ORDERING PHYSICIAN Not Given  
MEDICAL FACILITY Not Given  
ADDITIONAL RECIPIENT Not Given  
MEDICAL FACILITY ID Not Given  
PATHOLOGIST Not Given

**SPECIMEN**

SPECIMEN SITE Not Given  
SPECIMEN ID Not Given  
SPECIMEN TYPE Not Given  
DATE OF COLLECTION Not Given  
SPECIMEN RECEIVED Not Given

**Biomarker Findings**  
**Microsatellite status - MS-Stable**  
**Tumor Mutation Burden - TMB-Intermediate (11 Muts/Mb)**

**Genomic Findings**  
*For a complete list of the genes assayed, please refer to the appendix.*  
**EGFR amplification, L858R**  
**PTCH1 T416S**  
**CDKN2A/B loss**  
**RBMT0 Q494\***  
**TP53 R267P**

6 Disease relevant genes with no reportable alterations : **KRAS, ALK, BRAF, MET, RET, ERBB2, ROS1** 1

7 Therapies with Clinical Benefit in patient's tumor type    18 Clinical Trials  
 7 Therapies with Clinical Benefit in other tumor type

**BIOMARKER FINDINGS**

**Tumor Mutation Burden - TMB-Intermediate (11 Muts/Mb)**

3
 9 Trials see p.14

**Microsatellite status - MS-Stable**

**GENOMIC FINDINGS**

**EGFR - amplification, L858R**

3
 4 Trials see p.15

**PTCH1 - T416S**

THERAPIES WITH CLINICAL BENEFIT (IN PATIENT'S TUMOR TYPE)	THERAPIES WITH CLINICAL BENEFIT (IN OTHER TUMOR TYPE)
Atezolizumab	Avelumab
Nivolumab	Durvalumab
Pembrolizumab	

No therapies or clinical trials. see Biomarker Findings section

THERAPIES WITH CLINICAL BENEFIT (IN PATIENT'S TUMOR TYPE)	THERAPIES WITH CLINICAL BENEFIT (IN OTHER TUMOR TYPE)
Afatinib	Cetuximab
Erlotinib	Lapatinib
Gefitinib	Panitumumab
Osimertinib	
none	Sonidegib

- 1
**Pertinent Negatives**  
 Identifies important negative results that can be used for patient management
- 2
**Therapies with Clinical Benefit**  
 Interpretive content that can be used for patient management according to professional guidelines in oncology
- 3
**Clinical Trials**  
 Identifies trials based on patients' unique genomic profile with page number for quick reference

<b>TO LEARN MORE:</b>	<b>TO ORDER:</b>
Visit <a href="http://www.foundationmedicine.com/f1cdx">www.foundationmedicine.com/f1cdx</a>	Create an account to order online at <a href="http://www.foundationmedicine.com/signup">www.foundationmedicine.com/signup</a>

FoundationOne CDx™ is a next-generation sequencing based *in vitro* diagnostic device for detection of substitutions, insertion and deletion alterations, and copy number alterations in 324 genes and select gene rearrangements, as well as genomic signatures including microsatellite instability (MSI) and tumor mutational burden (TMB) using DNA isolated from formalin-fixed, paraffin-embedded (FFPE) tumor tissue specimens. For the complete intended use statement, including companion diagnostic indications, please see the FoundationOne CDx Technical Information, [www.foundationmedicine.com/f1cdx](http://www.foundationmedicine.com/f1cdx).